











How to Use a Pediatric Urine Collection Bag

For Girls	For Boys
<p>Step 1</p>  <p>Clean from front to back, cleaning between each skin fold of the labia.</p> <p>Clean around the anus.</p> <p>Allow area to air-dry.</p>	<p>Step 1</p>  <p>Start with cleaning the penis from the tip downward (pull back the foreskin if needed).</p> <p>Clean the scrotum. Clean around the anus. Allow area to air-dry.</p>
<p>Step 2</p>  <p>Remove the protective backing from the bottom half of the adhesive patch of the urine bag.</p>	<p>Step 2</p>  <p>Remove the protective backing from the bottom half of the adhesive patch of the urine bag.</p> <p>If your boy is active, it might be easier to leave the adhesive backing in place until you have fitted the opening over his scrotum and penis.</p>
<p>Step 3</p>  <p>Spread your child's hips and stretch out the skin folds in the groin.</p> <p>Put the bottom half of the bag on your child. Start with placing the middle on the tiny area between the anus and vagina. Work outwards, pressing the adhesive patch smoothly and firmly in place.</p>	<p>Step 3</p>  <p>Spread your child's hips and stretch out the skin folds in the groin.</p> <p>Put the bottom half of the bag on your child. Start with placing the middle on the tiny area between the anus and scrotum. Work outwards, pressing the adhesive patch smoothly and firmly in place around the scrotum.</p> <p>Place the penis into the opening of the urine bag.</p>
<p>Step 4</p>  <p>Remove the upper half of the adhesive patch.</p> <p>Working upward and outward, press the adhesive patch smoothly and firmly in place.</p> 	<p>Step 4</p>  <p>Remove the upper half of the adhesive patch.</p> <p>Working upward and outward, press the adhesive patch smoothly and firmly in place.</p> 



PEIRSON CENTER FOR CHILDREN

9820 E Burnside Street, Portland, OR 97220

971.202.7356